







CARACTERÍSTICAS
FEATURES

GROSOR

TOTAL VISCO 2CM

HR35 KG

VISCO 50KG/M3
THERM-ESSENCE

CUALQUIER BASE
ALL TYPES BASES

ERGONÓMICO
ERGONOMIC

7 ZONAS DESCANSO
7 REST ZONES

TRANSPIRABLE
BRAETHABLE

ORTOPÉDICO
ORTHOPEDIC

FUNDA IMPERMEABLE

FUNDA COVER

ANTIÁCAROS

ECOLÓGICO


CÁLIDO
WARM

FRESCO
FRESH

SUAVE
SOFT

FIRME
FIRM

COMPOSICIÓN
COMPOSITION



1 Tejido (Punto PES % 105gr/m2 lamina 70gr/m2 PU (Total 175gr/m2))	Fabric
2 Visco Therm 50KG 2cm	Memory Therm 50kg 2cm
3 HR 35kg 13cm	HR 35 kg 13 cm
4 Tejido (Punto PES % 105gr/m2 lamina 70gr/m2 PU (Total 175gr/m2))	Fabric

INFORMACIÓN TÉCNICA

- Impermeable / Transpirable
- Ignífugo crib 5 (BS7175)
- Acabado antibacteria + antihongos
- Lavable 90° C no planchar
- Pegado cola base agua (libre de disolventes)
- Funda con solapa

NUCLEO PERFILADO / GERIATRICO

